


jGCaMP8: Ultra-sensitive protein calcium sensors with fast kinetics

The [Looger Lab](#) and the [GENIE Project Team](#) at HHMI Janelia have developed a new suite of jGCaMP8 Calcium indicators, built upon the GCaMP scaffold. The jGCaMP8 sensors have fast kinetics without compromising sensitivity, setting a new standard for *in vivo* imaging. Sensors that have been extensively tested in mammalian neurons *in vivo* and *in vitro* are:


- **jGCaMP8f** (*fast*): 4x faster rise time, 2.5x faster decay time than jGCaMP7f
- **jGCaMP8m** (*medium*): almost 4x faster rise time and 3.5x more sensitive than jGCaMP7f
- **jGCaMP8s** (*sensitive*): 2x more sensitive than jGCaMP7s, >2x faster than jGCaMP7f (at 1 AP)

In vivo performance


Below are representative traces based on simultaneous fluorescence and cell-attached electrophysiology in mouse visual cortex. Small fields of view were imaged *in vivo* at effective time resolutions of 2 ms



Grand average of calcium transients elicited by single action potentials *in vivo*, aligned to action potential peak (red bar), resampled at 500 Hz. Right: zoom-in to rise kinetics



Cell-wise analysis of calcium transients elicited by single action potentials *in vivo*. Black dots represent single cells, bars represent median of single cells, horizontal lines represent values calculated from the grand averages (calculated from previous figure).



(Previous page) Representative *in vivo* ground truth recordings of jGCaMP8 sensors. Fluorescence traces were filtered with a Gaussian filter ($\sigma = 5$ ms).

Cultured neuron screen performance

The following table summarizes the sensitivity and kinetics of the jGCaMP8 sensors in the GENIE field stimulation-based [cultured neuron screen](#).

	Max dF/F	Half-rise time (ms)	Time to peak (ms)	Half-decay time (ms)
jGCaMP7f (control)	0.21±0.1	24.8±6.6	99.5±30.2	181.9±76.0
jGCaMP8f	0.41±0.12	7.1±0.74	24.8±6.1	67.4±11.2
jGCaMP8m	0.76±0.22	7.1±0.61	29.0±11.2	118.3±13.2
jGCaMP8s	1.11±0.22	10.1±0.86	57.0±12.9	306.7±32.2
jGCaMP8.712*	0.66±0.18	10.9±1.24	41.6±8.1	94.8±13.3

*jGCaMP8.712 is currently being tested and may be released at a later date

Characterization in purified protein

The jGCaMP8 indicators were purified for *in vitro* characterization of affinity, kinetics, and absorption.

	Kd (nM)	Hill coeff.	Max dF/F	k _{off} (fast)	k _{off} (fast) %	k _{off} (slow)	k _{off} (slow) %
jGCaMP7f (control)	150±2	3.10±0.16	31.0±1.1	7.34±0.12			
jGCaMP8f	334±18	2.08±0.22	78.8±9.7	37.03±0.75	91	1.37±3.49	9
jGCaMP8m	108±3	1.92±0.12	45.7±0.9	18.25±0.9			
jGCaMP8s	46±1	2.20±0.13	49.5±0.1	3.68±0.1			
jGCaMP8.712*	292±16	2.17±0.15	63.6±4.1	18.8±0.44			

*jGCaMP8.712 is currently being tested and may be released at a later date

	pKa, apo	pKa, sat	EC _{sat} (M ⁻¹ cm ⁻¹)	EC _{apo} (M ⁻¹ cm ⁻¹)
jGCaMP7f (control)	8.68±0.13	6.71±0.06	52458	2783
jGCaMP8f	7.71±0.06	6.68±0.01	50759	1927
jGCaMP8m	7.40±0.02	6.68±0.07	49856	2249
jGCaMP8s	7.65±0.04	6.51±0.04	56960	2116
jGCaMP8.712*	9.01±0.13	6.83±0.06	41862	1358

*jGCaMP8.712 is currently being tested and may be released at a later date

[Order jGCaMP8 on Addgene](#)

[Visit our website for more information](#)

[Contact us](#)

If you use these reagents in your work, please cite this document:

Y. Zhang, M. Rózsa, D. Bushey, J. Zheng, D. Reep, G. J. Broussard, A. Tsang, G. Tsegaye, R. Patel, S. Narayan, J. X. Lim, R. Zhang, M. B. Ahrens, G. C. Turner, S. S.-H. Wang, K. Svoboda, W. Korff, E. R. Schreiter, J. P. Hasseman, I. Kolb, L. L. Looger, *jGCaMP8 Fast Genetically Encoded Calcium Indicators* (2020). doi:10.25378/janelia.13148243.